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

www.ikts.fraunhofer.de

FRAUNHOFER INSTITUTE FOR CERAMIC TECHNOLOGIES AND SYSTEMS IKTS

- 1 Two-channel sensor platform for refractive index measurements; implementation with microfluidic system.
- 2 Electron microscope image of the nanostructured metallic sensor surface.

PROCESS-INTEGRATED REFRACTIVE INDEX SENSOR

Refractive index measurement in liquids

The refractive index is a standardized parameter in different production processes in chemistry and biotechnology. It allows for instance to evaluate the concentration of certain components in a liquid. However, in-line measurements with process-integrated refractive index sensors are still rarely seen.

Plasmonic sensor effect

The system analyzes the spectral position of the surface plasmon resonance in a nanostructured metal surface which is sensitive to the refractive index it is in contact with. For multisensor purposes the refractive index sensor array can be flexibly equipped with different surface functionalization routines. This allows for example the measurement of biomolecular concentrations by means of antibody-based surface modifications which induce a refractive index shift.

A robust sensor platform

The refractive index sensor platform developed at Fraunhofer IKTS features a high degree of miniaturization, a robust design, and the ability for realtime and multiplexed measurements. It is therefore well suited for reliable process integration – either via a microfluidic system or directly as a probe in contact with the liquid of interest. The current two-channel system has dimensions of $7 \times 7 \times 7\text{cm}^3$ and provides a resolution of the refractive index of 10^{-4} .

The system is set up in a modular way and thus allows addressing customer-specific adaption.

In addition, test and service measurements can be performed at Fraunhofer IKTS.

Fraunhofer Institute for Ceramic Technologies and Systems IKTS

Maria-Reiche-Strasse 2
01109 Dresden, Germany

Contact

Dr. Thomas Härtling
Phone +49 351 88815-550
thomas.haertling@ikts.fraunhofer.de

www.ikts.fraunhofer.de